Titration Master - User Manual & Tutorial

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Welcome to **Titration Master**, a virtual chemistry laboratory designed to explore the fascinating world of acid-base titrations. This interactive simulation allows you to perform experiments with various acids, monitor pH changes in real-time, and master the art of chemical analysis.

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Main Chart / Gráfica principal

- Displays the **Titration Curve** (pH vs. V_{base}).
- Highlights **Equivalence Points** (dots for polyprotic acids).
 n
- Shows the **Transition Range** of the selected indicator as a shaded band.

Live View / Vista en vivo

- Buret & Beaker:** Visual representation of the experiment.
- Liquid Color:** Changes dynamically based on the current pH and selected indicator.
- Digital Displays:** Real-time pH, volume added, and substance identification.
- Controls:** Use the progress slider or the **+DROP** button for precision.

2. Laboratory Modes / Modos de Laboratorio

Free Exploration / Exploración Libre

Experiment with any acid and concentration. Perfect for visualizing the difference between strong acids (HCl) and weak acids (CH_3COOH), or observing the multiple “steps” in polyprotic acids like Phosphoric Acid (H_3PO_4).

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The identity of the acid is hidden. You must perform the titration, analyze the shape of the curve and the equivalence points, and identify the unknown substance.

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You are given a **Target pH**. Your goal is to adjust the titration progress until the solution reaches exactly that pH (within ± 0.02 units).

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3. Step-by-Step Tutorial / Titorial paso a paso

Preparation / Preparación

1. Set the **Mode** to “Free Exploration”.
2. Select **Acetic Acid** (CH_3COOH).
3. Set $C_a = 0.10$ M and $V_a = 25$ mL.
4. Choose **Phenolphthalein** as the indicator.

The Titration / A Valoración

1. Slowly move the **Titration Progress** slider.
2. Watch the pH rise slowly as you add NaOH.
3. As you approach V_m (when $V_a = V_m$), the pH will jump sharply.
4. Observe the beaker: at $pH \approx 8.2$, the liquid starts turning pink.
5. Use the **+DROP** button for fine control near the equivalence point.

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pH and K_a

The value pK_a indicates the strength of an acid. A lower pK_a means a stronger acid.

- **Strong Acids** (HCl, H_2SO_4): Fully dissociate. The equivalence point with NaOH is at $pH = 7$.
- **Weak Acids**: Partially dissociate. The equivalence point is usually at $pH > 7$ due to the basic nature of the conjugate base.

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Polyprotic Acids / Ácidos Polipróticos

Acids like **Sulfuric** (H_2SO_4), **Oxalic** ($H_2C_2O_4$), or **Phosphoric** (H_3PO_4) have multiple acidic protons. Each proton has its own pK_a and its own equivalence point on the curve.

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5. Tips for Success / Consellos para o éxito

Indicator Choice: For strong acid/strong base titrations, almost any indicator works. For weak acids, you need an indicator like Phenolphthalein that changes in the basic range. **Precision:** Near the equivalence point, the pH changes very rapidly with very small volumes of base. This is the most critical part of the experiment!

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